

JUMPER

# FINGERTIP PULSE OXIMETER

Product model: JPD-500C (OLED)



## FEATURES

JUMPER



SpO2 Measuring



Pulse Rate Measuring



Adjustable Brightness



Multiple Display Modes




Alert Function



Automatic Power-Off

## STRUCTURE

JUMPER



## PACKAGING

JUMPER

CTN SIZE(CM)	QTY/CTN	GW(KG)	FCL/20
41 x 38 x 27	100	10	600

### ACCESSORIES

Hanging rope × 1  
User manual × 1

Dry batteries × 2



## SPECIFICATION

JUMPER

Display	OLED Screen
Direction	Four directions display
Parameters	Pulse Rate, SpO2
Accuracy	SpO2: 70%-100%, ± 2% Pulse Rate: 25-250bpm, ± 2bpm
Indication	Real-time Battery Status, Value Limited Remind
Application Range	Adult, Adolescent
Battery	AAA × 2
Dimensions	58.0 × 32.0 × 33.2 mm
Net Weight	28.6 g