

INJECTION TRAINER

00310 Simulated Patient



Our soft tissue injection pad is designed for practising intradermal, subcutaneous and intramuscular tissue injection techniques.

The Injection Trainer has multiple tissue layers representing the epidermis, dermis, fat and muscle layer, and can easily attach to an arm or thigh to help teach professional-to-patient communication.

SKILLS

- Subcutaneous injection
- Intradermal and intramuscular injection
- Management of tissue
- Professional-to-patient communication

QUALITIES

ANATOMY

– Epidermis, dermis, fat and muscle

REALISM

– Tissues are soft and warm to the touch

KEY FEATURES

– Epidermis layer peels back to quickly release subcuticular liquid

VERSATILITY

– Strap for hybrid simulation

SAFETY

– Latex free

CLEANING & CARE

– Durable replaceable epidermis

CONTAINS

- **00311** Skin Pad & Muscle Block
- **00320** Epidermis for Injection Trainer (x2)
- Strap-on base plate
- Soft tissue tray with liner