

KOKEN Breastfeeding Simulation Set

LM-113A



Outline

The optimal set for experiencing breastfeeding support.

Users can feel adequate close fit when wearing the KOKEN Postpartum Breast Care Model and can take a natural breastfeeding position.

By attaching the Latch on mask for KOKEN Baby to their KOKEN Baby, users can simulate a baby's appropriate latch.



KOKEN Postpartum Breast Care Model LM-113



Latch on mask for KOKEN Baby LM-113B



* KOKEN Baby is not included

Features

- By using both the KOKEN Breastfeeding Simulation Set and KOKEN Baby together, users can practice a "handsoff" approach to the support of breastfeeding.
- Users can observe and practice caring for the nipples and breasts.
- Simulation of engorgement of the breast enables users to observe the engorged state and practice methods for alleviating engorgement.
- Users can learn through experience that a correct latch cannot be achieved when the breast is engorged.



Components

KOKEN Postpartum Breast Care Model

Harness	1
Air pump	1
Baby powder	1
Protective board	1
Breast bag (right and left)	1

Latch on mask for KOKEN Baby

Latch on mask for KOKEN Baby	1
------------------------------	---

Specifications

KOKEN Postpartum Breast Care Model

Weight	Approx. 1.6kg
--------	---------------

KOKEN Postpartum Breast Care Model

Size	Apporox. 15(L) × 13(W) × 12(H)cm
------	----------------------------------

Weight	Approx. 96g
--------	-------------