



Prestan® Professional Adult Manikin

Adult Manikin features:

- Realistic looking and performing manikin is more acceptable to both instructors and students
- Far easier and faster to carry and setup at training sites
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training
- Face shield is secured so it won't slip and there's no need for adhesives
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway
- Located in the shoulder is an optional CPR Rate Monitor that provides instant visual feedback on compression rates - students know they are at the correct compression rate when two green lights illuminate

- The unique *Clicker Mechanism* provides an audible confirmation when the correct compression depth is reached, guiding to the appropriate force when administering CPR compressions.
- Available individually or in convenient multipacks of four
- Prestan Manikins are latex free



ITEM # PP-AM-400M-MS (with CPR Monitor)
ITEM # PP-AM-400-MS (without CPR Monitor)
(Prestan Professional Adult Manikin 4-Pack)

Reorder Accessories:

- #PP-ALB-10
(10-pack Adult Manikin face-shield/lung-bags)
- #PP-ALB-50
(50-pack Adult Manikin face-shield/lung-bags)
- #PP-AFS-50
(50-pack Adult Manikin face-shields)



ITEM # PP-AM-100M-MS (with CPR Monitor)
ITEM # PP-AM-100-MS (without CPR Monitor)
(Prestan Professional Adult Manikin Single)