

ULTRASOUND VASCULAR ACCESS TRAINER

60530



An echolucent and contoured slab with representation of the jugular vein for teaching needle insertion under ultrasound.

SKILLS

- Use of ultrasound technology
- Handling ultrasound transducers
- Insertion of needles under ultrasound guidance
- Handling and management of guidelines and catheters under ultrasound
- Recognition of vessels using ultrasound

QUALITIES

REALISM

- Echolucent material provides realistic ultrasound image

KEY FEATURES

- Vein is self-sealing for repeated use and rechargeable with water via a one-way valve
- Designed to be used on the bench top

VERSATILITY

- Fluid volume and pressure can be changed in the vessel to represent variability in patient circulation, providing various levels of difficulty for needle insertion

SAFETY

- Veins contain latex

CONTAINS

- **60531** Ultrasound Vascular Access Slab
- **60532** Ultrasound Vascular Access Vein
- Base
- 10ml syringe