

Oxygen Concentrator



8F-3W/8F-5W Oxygen Concentrator



Global sales leader



Medical grade



95.5%-87% oxygen concentration



ultra-silence 41dB

Flow Range: 0.5-5L/min

Oxygen concentration: 93% ± 3%

6.5 inch screen with remote

Maximum Outlet Pressure: 40-70kPa

Sound Level: ≤39dB(A)

Weight: ≤25kg

Dimension: 48cm × 30cm × 65cm



8F-5AW Oxygen Concentrator



Global sales leader



Medical grade



95.5%-87% oxygen concentration



5L oxygen flow

Flow Range: 0.5 ~ 5L/min

Oxygen Concentration: 95.5%-87%

Maximum Outlet Pressure: 40 ~ 70kPa

Sound Level: Front side 48dB(A)(average), Whole unit 52dB(A)

Input Power: 400VA

Weight: 15.5kg (34.1 lbs)

Dimension: H50cm(19.7") × W39cm(15.4") × D 24.5cm(9.6")

Minimum Operating Time: 30minutes



9F-5AW Oxygen Concentrator

Flow range: 0.5 ~ 5L/min

Oxygen concentration: 93% ± 3%

Output pressure: 40~70kPa

Input power: 400VA

N.W.: 18kg

Sound level: Front side 49.5 dB(A), whole unit 52 dB(A)

Size: 50*21*64cm





9F-3AW Oxygen Concentrator

Flow range: 0.5 ~ 3L/min
 Oxygen concentration: 93% ± 3%
 Output pressure: 20~50kPa
 Input power: 320VA
 N.W.: 16.5kg
 Operation noise level: ≤ 46dB(A)
 Size: 50*21*64cm

- High intake position, cleaner air
- Remote control function
- 4 alarm system (Low oxygen concentrator alarm, Power outage alarm, Pressure cycle failure alarm, Compressor failure alarm) , self detection of fault
- 4 safety protection

Oil free pump

-The heart of Oxygen Concentrator

- Compressor adopts with Germany technology
- Over 10000 hours service life
- Germany Technology



9F-3BW Oxygen Concentrator

-  Stable Oxygen Concentration
-  Backlight & brightness adjustable LCD screen
-  Ultra Silent
-  Infrared Remote Control System
-  Voice broadcast function
-  U.S. IMPORTED SIEVE

Flow range: 0.5 ~ 3L/min
 Oxygen concentration: 93% ± 3%
 Output pressure: 20 ~ 50kPa
 Input power: 230VA
 N.W.: 10.5kg
 Operation noise level: ≤ 46dB(A)
 Size: 31.5*20*39.5cm



YU600 Homecare Oxygen Concentrator

Flow range: 1-3L/min
 Oxygen concentration: 70%-90%
 Input power: 150VA
 N.W.: 8kg
 Operation noise level: 45dB(A)
 Size: 31.7*19*37.7(cm)



YU300 / YU500 Homecare Oxygen Concentrator

Flow range: 1 ~ 5L/min
 Oxygen concentration: 30%~90%
 Input power: 100VA
 N.W.: 7.5kg
 Operation noise level: 43dB(A)
 Size: 27cm*17.5cm*29.5cm